

uBase™

iliosacral fixation system

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100 years of reliability and German precision

Spinal systems by Ulrich Medical® stand for quality “Made in Germany.”

They are the successful result of systematic development activities and many years of experience in medical technology. For more than 100 years, we have done our best every day for our customers and for patient health with our proprietary and innovative products.

Leading-edge technology and competence from one source

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General information

Please ensure that you are using and observing the most current instructions for use and surgical technique of this system at all times. Please note also the “Instructions for use” WS 714695 T-handle with torque limitation, 9 Nm, coupling V (CS 7146). These can be downloaded any time free of charge at: www.ifu.ulrichmedical.com.

The information regarding sterile implants are valid exclusively for countries in which the implants are approved. All unsterilized implants are also available in sterile packaging. To order sterile products, add “-S” to the article number.



This surgical technique describes the implant, the instruments, and the steps involved in the application of the uBase™ system in the case of open access. It is not sufficient as the sole basis for a successful application of the system. It is recommended to master the surgical technique with an experienced surgeon. Please note the “Instructions for use” for uBase™ implants.

Each uBase™ screw is cannulated and thus permits secure and precise placement via guide wires. The screws can be implanted in the sacrum as well as in the ilium and also in connection with Brainlab navigation if needed.

The uBase™ system is characterized by a low-profile screw head which enables a stable-angle connection to the rod using connectors provided. This reduces the need for rod contouring and yields mechanically stable instrumentations.

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The small pictures at the bottom of the page show the step-by-step application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-colored.

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Intended use

uBase™ is an implant system for the posterior surgical stabilization and fixation of the human sacrum and ilium in combination with uCentum™ or tangoRS™.

Indications

Spinal instrumentations which require additional stabilization and fixation of the sacrum and/or ilium, taking the indications and contraindications of the fixator used into account.

Contraindications

- Patients with fever or leukocytosis in the case of non-spine associated infections
- Patients with a documented history of material allergy or who tend to react to foreign bodies

**MR conditional**

Non-clinical testing has demonstrated that the implants of the uBase™ system are MR conditional. For further information please refer to the respective IFU.

**Additional information for processing, assembly and disassembly**

Please follow UH 1100 "Processing manual implants and instruments". For individual instruments which are correspondingly marked in the surgical technique (⚠), we refer to the "Assembly and disassembly instructions with special cleaning instructions". You can download them any time free of charge at: www.ifu.ulrichmedical.com.

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
Note: All uBase™ iliosacral screws are cannulated. It is recommended to implant them using the guide wire.

**CS 7171-01**

Montagestift

Easy-fit

**CS 7171-02**

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

Inlay for connector

**CS 7171-03**

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Nut for easy-fit

**CS 7171-06-01**Stab-Verbinder,
Ø 6,0 mm, kleinConnector for rod,
Ø 6.0 mm, small**CS 7171-06-02**Stab-Verbinder,
Ø 6,0 mm, mittelConnector for rod,
Ø 6.0 mm, medium**CS 7171-06-03**Stab-Verbinder,
Ø 6,0 mm, großConnector for rod,
Ø 6.0 mm, large

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Threaded element for
guidewire holder

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T-handle with torque
limitation, 9 Nm
coupling V

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Note: Please note also
the “Instructions for use”
WS 714695 T-handle with
torque limitation, 9 Nm,
coupling V (CS 7146).
These can be downloaded
any time free of charge at:
www.ifu.ulrichmedical.com.

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Gearshift handle M8, for navigation, coupling A

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Awl, cannulated, for IS-screw, Ø 7.5 mm, for navigation

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Awl, cannulated, for IS-screw, Ø 8.5 mm, for navigation

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Awl, cannulated, for IS-screw, Ø 10.0 mm, for navigation

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Drill, cannulated, for
IS-screw, Ø 7.5 mm,
coupling AO

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Drill, cannulated, for
IS-screw, Ø 8.5 mm,
coupling AO

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Drill, cannulated, for
IS-screw, Ø 10.0 mm,
coupling AO

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Tap, for IS-screw, Ø 7.5 mm,
cannulated, coupling A

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Tap, for IS-screw, Ø 8.5 mm, cannulated, coupling A

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Tap, for IS-screw, Ø 10.0 mm, cannulated, coupling A

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Screwdriver for IS-screw, self-retaining

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Adjuster for screwdriver for IS-screw, self-retaining

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Locking rod for screwdriver for IS-screw, self-retaining

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Screwdriver for IS-screw, closed, coupling A

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Screwdriver for nut, self-retaining

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Ball tip probe,
length 270 mm

55830-25A*

Brainlab Instrument star unit
ML (calibration with ICM4)



55830-29*

Brainlab Instrument star unit
L (calibration with ICM4)

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

**Purely sacral
 screw placement
 (S2 alar screw)**

The point of entry is here between an imaginary connecting line from the posterior sacral foramina of S1 and S2 and the lateral sacral crest. The screw position is frequently oriented parallel to the sacroiliac joint and ends at the anterolateral tip of the S2 ala. The sacroiliac joint should not be touched. In the case of bicortical screw placement, the correct screw length should be selected to avoid touching the iliac vessels and the presacral nerves.

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


**Sacroiliac screw
 placement
 (S2 alar-iliac screw)**

Screw entry is 2 to 4 mm lateral and 4 to 8 mm distal to the posterior S1 sacral foramen. The screw position is oriented on the curvature of the pelvis and is directed towards the anterior inferior iliac spine.

Caution: It is urgently recommended to study the literature beforehand.

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
**Purely iliac
 screw placement
 (iliac screw)**

Screw entry is lateral to the sacroiliac joint or at the posterior superior iliac spine. The tip of the screw either lies above the hip joint (select screw length such that it is not touched!) or in the direction of the anterior inferior iliac spine.


Caution: In the case of bicortical screw placement, there is a risk of injury to the gluteal arteries and gluteal nerves; the screw length should therefore be measured correctly.

Note: The use of computer navigation is recommended for the three placement techniques described.


Important: Pedicle screws may only be implanted in S1 in a convergent fashion.

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
General instrumentation information

The special design of the stable-angle iliosacral screw offers the user a high degree of flexibility when connecting the rod with the screw and enables contouring of the rod to be minimized. The connector itself has polyaxial movement of $\pm 8^\circ$. In addition, the following adaptations can be made:

Distance between screw and rod

To adjust the distance between the screw and the rod, three connectors of different lengths are available [1].

Equalization of height offset

By rotating the screw 180° , the position of the connector shifts in a vertical direction [2].

Alignment of the connector in the sagittal plane

By rotating the screw, the connector can be tilted such that it fits the length of the rod in the sagittal plane [3].

Offset connector from the uCentum™ system

For the case in which all of the above adaptations are not successful in establishing a connection, offset connectors (CS 3812-01-015, -030, -060) are available [4]. The application can be found in the uCentum™ surgical technique.

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Handling the instrument set and the introduction of the implants are shown using a model of the spine.

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Exposure of the spine

The appropriate open approach is selected and the posterior spinal region is exposed.

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1. Preparation using a cannulated awl and trocar wire

Assembling the cannulated awls

For the various IS screw diameters 7.5 mm, 8.5 mm and 10 mm, there are appropriate cannulated awls available (CS 7174-075, -085, -100).


The gearshift handle (CS 7174-02) is positioned on the cannulated awl and the trocar wire for navigation, sharp (CS 7174-01) is screwed in.

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
Preparing the screw hole

The cortical bone of the vertebral body is opened at the entry point using the pre-assembled cannulated awl for IS screws. The screw hole is prepared further, whereby the awl is driven further in a clockwise direction using constant rotation until it reaches the later intended position of the tip of the screw. Lateral bending of the awl should be avoided.

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The trocar wire is subsequently unscrewed and removed [2] and the awl remains in the sacrum/ilium. The awl can likewise be used to determine screw length [1].

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
Inserting the guide wire

The guide wire (CS 2624) is introduced through the cannulation of the awl [1]. The awl is then carefully removed by rotating clockwise, while the guide wire remains in situ [2].

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
Palpating the drilling channel

Using the ball tip probe (UL 8522 -27), the integrity of the screw hole can then be checked. The guide wire must be temporarily removed to do this [3].

Tapping (optional)

As support for the preparation in the case of sclerotic cancellous bone, the tap (CS 7176-xxx) can be used. The length of the tapping area of the instrument allows preparation of the screw hole over the entire screw length [4].

Important: The guide wire must be held in position so that it is not pushed forward. Lateral C-arm control.

If desired, the guide wire holder can also be used to hold the guide wire in position during tapping. The application is explained on pages 31 and 32.

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2. Navigation of the cannulated awl with Brainlab®


Analogously to the preparation of the screw hole described above using a cannulated awl and trocar wire, the procedure can also be performed in combination with Brainlab® navigation.

Attaching the reference star for use with Brainlab®

The Brainlab® reference star unit (55830-25A or 55830-29) must be fully mounted on the adapter on the awl. It must be mounted in an oriented fashion so that the pin on the awl engages in the corresponding slot of the reference star and the self-retaining mechanism snaps into place in the corresponding groove of the awl. Then the gearshift handle M8, for navigation (CS 7174-02) can be positioned on the awl and the trocar wire for navigation, sharp (CS 7174-01) can be screwed in.

The application can be performed using two reference stars which are to be calibrated independently of each other in combination with Brainlab navigation software and navigated together with different instruments (e.g. awl and probe).

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Note: Because of the long cutting lengths of the size 8.5 and 10 awls, the associated calibration sleeves must be used for the calibration of these awls in combination with Brainlab ICM4.

Mounting the calibration sleeve on the awl

The calibration sleeve for Brainlab® ICM4 (CS 7174-085-01 and CS 7174-100-01) must be pushed onto the tip until the self-retaining mechanism snaps into place.

Preparing the navigation

After the Brainlab® marker spheres have been screwed on, the instrument is ready for manual calibration. Using the Brainlab® Instrument Calibration Matrix (ICM4), the instrument can be easily calibrated. The instrument with reference star and handle must be held during use such that the navigation camera has a direct and unimpeded view of the reference star.

The awls with the calibration sleeve must be calibrated using the Brainlab® Instrument Calibration Matrix (ICM4) and Brainlab® navigation software. To do this, the awl is inserted in the V groove of the ICM4 and advanced all the way to the tip.

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Note: Detailed information on handling the Brainlab® Disposable Reflective Marker Spheres, the ICM4, the reference stars, manual calibration and navigation of the instruments in combination with the navigation software can be found in the respective Brainlab user manual.

Warning: The navigated instruments are very precise and very sensitive medical instruments and must be used with extreme care. If the instrument should fall or be damaged, the calibration must be carefully checked or the instrument must be sent out for inspection. Improper handling can cause serious injuries to the patient.

Warning: Plan the OR system setup and the position of the reference star on instruments and patients prior to the surgery. The navigation camera must have a direct and unimpeded view of all reference stars. The navigation system should be set up as described by Brainlab.

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Warning: The navigated awls (CS 7174-075, -085, -100) may only be used in connection with the associated gearshift handle M8 (CS 7174-02) and the trocar wires (CS 7174-01 and CS 7174-03). ulrich medical® is not liable for the use of non-specified instruments.

Warning: The cannulated and sharp awls may only be calibrated via the V-notch of the Brainlab® ICM4 and the associated calibration sleeves are to be used for this purpose for the size 8.5 and 10.0 cannulated awls.


Warning: After changing an awl, the instrument must once again be calibrated with the Brainlab® ICM4.

Warning: Ensure that bending forces on the calibrated instruments are avoided since this deviation can affect accuracy.

Warning: Do not change the orientation of the awl during use since this may cause the awl to break within the bone.

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3. Preparation with a awl and blunt awl

Opening the cortical bone

The cortical bone of the vertebral body is opened at the point of entry using the awl (CS 7174-04) together with a handle (CS 7139 or CS 7174-02) [1].

Note: Other handles from the uCentum™ (CS 7138) and tangoRS™ (CS 3018) systems are available upon request.

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
Assembling the cannulated awls

For the various screw diameters 7.5 mm, 8.5 mm and 10.0 mm, there are appropriate cannulated awls available (CS 7174-075, -085, -100). The gearshift handle (CS 7174-02) is positioned on the cannulated awl and the trocar wire, blunt (CS 7174-03) is screwed in [2].

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

Preparing the screw hole

The screw hole is prepared further, whereby the awl is driven further in a clockwise direction using constant rotation until it reaches the later intended position of the tip of the screw. Lateral bending of the awl should be avoided. The trocar wire is subsequently unscrewed and removed [3] and the awl remains in the sacrum/ilium. The awl can likewise be used to determine screw length [4].

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
Inserting the guide wire

The guide wire (CS 2624) is introduced through the cannulation of the awl [1]. The awl is then carefully removed by rotating clockwise while the guide wire remains in situ.

Note: To insert screws longer than 100 mm, the guide wire with a length of 580 mm (CS 7173-580) must be used.

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Palpating the drilling channel

Using the ball tip probe (UL 8522-27), the integrity of the screw hole can then be checked [2].

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Tapping (optional)

As support for the preparation in the case of sclerotic cancellous bone, the tap (CS 7176-xxx) can be used. The length of the tapping area of the instrument allows preparation of the screw hole over the entire screw length [3].

Important: The guide wire must be held in position so that it is not pushed forward. Lateral C-arm control.

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
If desired, the guide wire holder can also be used to hold the guide wire in position during tapping. The application is explained on pages 29 and 30.

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The cortical bone of the vertebral body is opened at the point of entry using the awl (CS 7174-04) together with a handle (CS 7139 or CS 7174-02) [1]. In the case of sclerotic bone, the guide wire can be initially positioned and then the screw hole can be prepared using a cannulated drill.

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If the guide wire cannot be positioned, non-cannulated drills are also available (CS 7174-075, -085, -100)

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The corresponding adapter for drill that fits the drill is screwed into the drill guide. For the various screw diameters 7.5 mm, 8.5 mm and 10.0 mm, there are appropriate cannulated drills available (CS 7175-075-01, -085-01, -100-01) [2]. The drill is connected to the drilling machine using the AO coupling.

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

Preparing the screw hole

The screw hole is further prepared with the drill in which the drill guide is used to help maintain position and direction during drilling. Then the drill is withdrawn [1].

| Schrauben-Durchmesser Diameter screw | Bohrer-Durchmesser Diameter drill | Anzahl der Ringe Number of rings |
|--|-------------------------------------|------------------------------------|
| 7,5 mm 7.5 mm | 4,8 mm 4.8 mm | 1 |
| 8,5 mm 8.5 mm | 5,6 mm 5.6 mm | 2 |
| 10,0 mm 10.0 mm | 6,8 mm 6.8 mm | 3 |

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Palpating the drilling channel

Using the ball tip probe (UL 8522-27), the integrity of the screw hole can then be checked [2]. Then the guide wire (CS 2624) is inserted in the screw hole.

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
Note: The screw length can be determined as follows:

- Using the scale on the awls
- Using the scale on the drill
- By palpating the drill hole with the ball tip probe and reading the length on the scale on the storage tray

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Tapping (optional)

As support for the preparation in the case of sclerotic cancellous bone, the tap (CS 7176-xxx) can be used. The length of the tapping area of the instrument allows preparation of the screw hole up to the opposing cortical bone.

Important: The guide wire must be held in position so that it is not pushed forward. Lateral C-arm control.

If desired, the guide wire holder can also be used to hold the guide wire in position during tapping [1]. The application is explained on pages 31 and 32.

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1. Insertion of the IS screws using the self-retaining screwdriver

Note: All uBase™ iliosacral screws are cannulated. It is recommended to implant them using the guide wire.

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


Mounting the self-retaining screwdriver for IS screw


First the desired handle, handle XL, Tr 12x3 (CS 7139) [1] or T handle with ratchet (CS 7140) [2] is mounted on the screwdriver for IS screw, self-retaining (CS 7177-01-01). The thumb wheel for screwdriver for IS screws (CS 7177-01-02) is inserted in the recess in the screwdriver [3]. Then the locking rod is guided from the front through the opening of the screwdriver up to the thumb wheel, taking note of the laser markings on the screwdriver and locking rod [4] and it is completely screwed in as far as it will go [6] by rotating the thumb wheel (in the "open" direction) [5].

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

Mounting the IS screw on the self-retaining screwdriver

The screw head of the screw to be implanted is held via the mounted self-retaining screwdriver [1] and guided into the recess [2] such that the laser markings on the head and instrument are lined up [3]. Then the screw is secured by rotating the thumb wheel as far as it will go whereby the orange marking in the viewing window at the thumb wheel [4] displays the correct fixation of the screw.


Implanting the IS screws

The mounted uBase™ iliosacral screw is advanced using the mounted screwdriver via the guide wire up to the bone and implanted [5]. In doing so, it should be ensured that the screw is not initially screwed in entirely up to the bone so that the screwdriver can still be removed from the screw head.

Important: The guide wire must be held in position so that it is not pushed forward when the screw is screwed in. Lateral C-arm control is recommended.

To loosen the secure connection between the screwdriver and the screw head, the screw head is once again released by rotating on the thumb wheel ("open" direction) and the screwdriver can be removed.

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To completely screw in the IS screws, the screwdriver for IS screw, closed (CS 7177-02) is used. The laser marking in the front on both screwdrivers can be used to monitor the position of the easy-fit thread in the screw head. The thread in each case is on the opposite side.

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

2. Insertion of the IS screws using the closed screwdriver

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First the screwdriver for IS screw, closed (CS 7177-02) is combined with the desired handle, handle XL, Tr 12x3 (CS 7139) [1] or T handle with ratchet (CS 7140) [2]. The screw to be implanted is advanced via the guide wire to the bone [3]. Then the closed screwdriver is also pushed over the guide wire and placed on the screw head so that the laser markings on the head and screwdriver line up [4]. The easy-fit thread is then on the opposite side. At the end, the guide wire can be removed.

Important: The guide wire must be held in position so that it is not pushed forward when the screw is screwed in. For an alternative solution with the guide wire holder, see page 31. Lateral C-arm control is recommended.

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
3. Implanting the screw with the guide wire holder (optional)

The guide wire holder prevents the guide wire from shifting during screw insertion. When correctly used, it is an important safety feature of the system. The use of the guide wire holder is recommended particularly in the case of bicortical screw placement.

Note: For the insertion of screws longer than 100 mm, ensure that the guide wire with a length of 580 mm (CS 7173-580) is used.

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
Assembly of the guide wire holder


Important: The guide wire holder can only be used with the handle XL, Tr 12x3 (CS 7139). Use with the T-handle with ratchet (CS 7140) is not possible!

The threaded element (CS 3822-01) is screwed into the handle XL as far as it will go [1]. Then the handle for the guide wire holder (CS 3822-02) is screwed into the threaded element [2]; in doing so, the handle should not be screwed all the way into the threaded element [3].

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
Screwing in the screw with the guide wire holder (optional)

The uBase™ IS screw is implanted via the guide wire in place using the screwdriver and mounted guide wire holder.




As soon as the tip of the screw comes into contact with bone, the guide wire is connected to the threaded element (CS 3822-01) by screwing in the handle (CS 3822-02) further [1].

Important: While the IS screw is screwed in, the handle of the guide wire holder must be held since it should not simultaneously be rotated [2]. In this way, the guide wire is held in position while the IS screw is screwed in [3].

To release the tight connection between the screwdriver and the screw head, the guide wire must initially be unscrewed out of the handle by rotating the guide wire holder and then guided out back. Thereafter, the screw head can be released once again by rotating the thumb wheel ("open" direction) and the screwdriver can be removed. If necessary, the IS screw can be screwed in on the bone using the screwdriver for IS screw, closed (CS 7177-02 and CS 7140). At the end, the guide wire can be removed.

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4. Completing the instrumentation

Inserting the easy-fit

The easy-fit (CS 7171-01) must initially be inserted in the screwdriver for easy-fit (CS 1461) [1+2]. Thereafter, the easy-fit can be screwed into the screw head [3].

Important: The easy-fit must be fully screwed into the screw head in order to ensure the correct position of the nominal breaking point for the nut.

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
Alignment of the screw heads

The easy-fit is aligned by rotating the screw head using the screw adjuster (CS 7178-03) until the easy-fit is vertical to the surface of the bone [4].

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
Connecting the IS screw with the rod

To connect the IS screw with the rod of the screw-rod system, rod connectors of varying lengths for Ø 6.0 mm rods (CS 7171-06-XX) are used. These are initially slid onto the rod. The forceps for connector (CS 7179) facilitate the sliding if the rod is already contoured. The appropriate connector is now placed over the easy-fit with the rod and positioned on the screw head. If necessary, the screw head must be aligned once more.

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
Option: If a connection cannot be established with the existing rod connectors, there is still the option of achieving a connection with the rod using an offset connector (CS 3812-01-15, -30, -60) with a length of 15, 30 or 60 mm from the uCentum™ system range. For this purpose, the offset connector is inserted in the rod connector. The application of the offset connector can be found in the surgical technique, open, of the uCentum™ systems.

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

Inserting the inlay and nut

The inlay for connector (CS 7171-02) is guided over the easy-fit and inserted into the connector [1].

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
Ensure that the easy-fit is screwed into the screw head as far as it will go to ensure the correct position of the nominal breaking point to the nut.

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

Securing the instrumentation

To secure rod connectors and the IS screw on the rod, the nut for easy-fit (CS 7171-03) is screwed in a self-retaining fashion on the easy-fit using the screwdriver for nut [2-4].

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

Then, using the screwdriver for nut (CS 7178-02) together with the handle with torque limitation, 9 Nm (CS 7146) [1], the instrumentation first undergoes final tightening. After that, the easy-fit is separated by lifting the screwdriver above the nut and tilting the instrument [2-4].

Important: After the separation procedure, the nut must once again be tightened using the handle for torque limitation, 9 Nm (CS 7146).

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
The uBase™ instrumentation is removed in the sequence opposite to the sequence of implantation steps.

First the T-handle with torque limitation, 9 Nm (CS 7146) is mounted using the screwdriver for nut (CS 7178-02). Then the torque screwdriver is used to remove the nut. The remaining part of the easy-fit must subsequently be removed using forceps. Then the rod connector can be removed from the screw head. The rod connector can be spread apart with the forceps for connector (CS 7179) so that it can be pushed from the rod more easily.


The screw is then removed using the closed screwdriver for IS screw (CS 7177-02 with CS 7140).

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
Layer 1 for uBase™
implants, for screws
Ø 7.5 mm

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
Layer 2 for uBase™
implants, for screws
Ø 10.0 mm

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Layer 3 for uBase™
implants, for screws
Ø 8.5 mm

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